



RETRACTOR MAKER

PRODUCT BROCHURE OF

Knapp Lacrimal Sac Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Knapp Lacrimal Sac Retractor

SKU RM-KLSR-1

\$39 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	1 Prong (Single Hook)	Available Widths: 5.0 mm · Available Shafts & Tips: Rigid or Flexible (Sharp / Blunt) · Price (USD): \$39.00 – \$49.00
1	2 Prongs	Available Widths: 5.0 mm or 6.0 mm · Available Shafts & Tips: Rigid or Flexible (Sharp / Blunt) · Price (USD): \$39.00 – \$49.00
1	3 Prongs	Available Widths: 6.0 mm or 8.0 mm · Available Shafts & Tips: Rigid or Flexible (Sharp / Blunt) · Price (USD): \$39.00 – \$49.00
1	4 Prongs	Available Widths: 6.0 mm or 8.0 mm · Available Shafts & Tips: Rigid or Flexible (Sharp / Blunt) · Price (USD): \$39.00 – \$49.00
1	Protective Sleeve / Box	Available Widths: Standard Transit Packaging · Available Shafts & Tips: — · Price (USD): Included

PRODUCT DESCRIPTION

The Knapp Lacrimal Sac Retractor is a specialized ophthalmic and oculoplastic surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement during dacryocystorhinostomy (DCR) and lacrimal apparatus procedures. Designed with an ultra-fine multi-prong rake head, this instrument allows surgeons to retract delicate medial canthal tissues cleanly, providing clear exposure of the lacrimal fossa while safeguarding the lacrimal sac from puncture.

Built with exceptional customizability, this Knapp Lacrimal Sac Retractor is offered in three overall lengths: 14 cm (\$39), 15.3 cm (\$44), and 16 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select between a Rigid Shaft for high stability during delicate maneuvers or a Flexible Shaft built to grant adaptive intraoperative maneuverability around complex facial angles. The rake head is configurable in 1, 2, 3, or 4 prongs with Sharp or Blunt tips across 5.0 mm, 6.0 mm, and 8.0 mm working widths.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Knapp Lacrimal Sac Retractor to strict clinical tolerances, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for oculoplastic, ENT, and microsurgical suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)