



RETRACTOR MAKER

PRODUCT BROCHURE OF

Fan Retractor With 5 Fingers

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Fan Retractor With 5 Fingers

SKU RM-FR5F-299-1

\$299 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
5 mm — Articulating	RM-FR5F-5A-1	\$299
5 mm — Non-Articulating	RM-FR5F-5N-1	\$299
10 mm — Articulating	RM-FR5F-10A-1	\$299
10 mm — Non-Articulating	RM-FR5F-10N-1	\$299

PRODUCT DESCRIPTION

The Fan Retractor With 5 Fingers is a premier laparoscopic surgical retractor engineered specifically to provide expansive, atraumatic organ displacement and broad-surface visceral elevation during minimally invasive general, bariatric, and gynecological procedures. Designed with five interlocking mechanical blades that collapse into a slender profile for trocar insertion, this instrument expands inside the peritoneal cavity into a broad fan-shaped lattice that holds large anatomical structures smoothly away from operative dissection planes.

Built to clinical framework specifications, this Fan Retractor With 5 Fingers model (\$299) is available in 5 mm or 10 mm shaft diameter options, with working lengths ranging from 330 mm to 375 mm (overall instrument length 438 mm to 487 mm). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select between Articulating or Non-Articulating shaft functional types to navigate around complex pelvic and abdominal obstacles. The cannulated chassis integrates a dedicated luer-lock flush port that ensures thorough lumen cleansing and biological residue removal. The smooth, rounded finger margins disperse counter-traction evenly to safeguard fragile visceral parenchyma.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Fan Retractor With 5 Fingers to withstand repetitive hospital autoclave cycles, providing laparoscopic operating theaters, bariatric suites, and digestive disease centers worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)