



RETRACTOR MAKER

PRODUCT BROCHURE OF

Senn Miller Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Senn Miller Retractor

SKU RM-SMR-SERIES-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
15.0 cm Reach — 10×9 mm flat blade / 10×11 mm 3-prong rake	RM-SMR-SERIES-1-1	\$29
15.5 cm Reach — 10×9 mm flat blade / 10×11 mm 3-prong rake	RM-SMR-SERIES-1-2	\$34
16.0 cm Reach — 10×9 mm flat blade / 10×11 mm 3-prong rake	RM-SMR-SERIES-1-3	\$39
20.0 cm Reach — 10×9 mm flat blade / 10×11 mm 3-prong rake	RM-SMR-SERIES-1-4	\$49

PRODUCT DESCRIPTION

The Senn Miller Retractor is an indispensable double-ended surgical retractor recognized worldwide for providing versatile, multi-depth soft-tissue displacement and wound margin exposure during general, plastic, and orthopedic procedures. Designed with an L-shaped flat blade on one end and a three-prong rake on the opposing end, this instrument allows surgical teams to transition effortlessly between broad tissue holding and non-slip skin traction without exchanging instruments.

Built to withstand frequent autoclave sterilization cycles while providing a lifetime warranty against rust, this Senn Miller Retractor series is offered across four overall lengths: 15 cm (\$29), 15.5 cm (\$34), 16 cm (\$39), and 20 cm (\$49). Handcrafted from premium operating room-grade materials, its integrated double-ended blade dimensions feature a flat/smooth blade end measuring 10 mm × 9 mm paired with a 3-prong rake end measuring 10 mm × 11 mm. Surgeons can customize the rake end with Blunt End Prongs engineered for safe, atraumatic retraction of soft tissues, or Sharp End Prongs engineered to deliver a firm, slippage-free grip on tougher tissue edges.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Senn Miller Retractor to uphold strict clinical tolerances, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for hospital operating rooms, outpatient clinics, and plastic surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)