



RETRACTOR MAKER

PRODUCT BROCHURE OF

McBurney Thyroid Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



McBurney Thyroid Retractor

SKU RM-MTR-19-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard — 7.5 in, 16x10 mm blade	RM-MTR-19-1-1	\$49

PRODUCT DESCRIPTION

The McBurney Thyroid Retractor is a specialized endocrine and cervical surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and strap muscle retraction during thyroid and parathyroid procedures. Designed with an anatomically curved right-angled blade, this instrument allows surgical teams to mobilize the sternohyoid and sternothyroid muscles smoothly away from the thyroid gland, establishing an unobstructed operative window while safeguarding adjacent neurovascular structures.

Built to ensure balanced mechanical leverage during delicate neck dissections, this dedicated McBurney Thyroid Retractor configuration features an overall tool length of 19 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to guarantee lasting corrosion resistance, its working end incorporates a 16 mm long by 10 mm wide blade profile with rounded, atraumatic margins that prevent localized tissue contusions or tears. The ergonomic flat handle provides operating assistants with steady manual leverage throughout extended operative cases.

Forged with an anti-glare satin matte finish, Retractor Maker produces the McBurney Thyroid Retractor to withstand repetitive hospital autoclave cycles, providing endocrine surgical suites, ENT operating theaters, and general surgical clinics worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)