



RETRACTOR MAKER

PRODUCT BROCHURE OF

Harrington Sweetheart Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Harrington Sweetheart Retractor

SKU RM-HSR-1

\$129 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
4 × 12.8 cm	RM-HSR-1-1	\$129
6.5 × 12.8 cm	RM-HSR-1-2	\$134
6.5 × 14.6 cm	RM-HSR-1-3	\$139
9 × 25.5 cm	RM-HSR-1-4	\$144
11.5 × 25.5 cm	RM-HSR-1-5	\$149

PRODUCT DESCRIPTION

The Harrington Sweetheart Retractor is a premier visceral surgical retractor engineered specifically to provide broad, atraumatic displacement of parenchymal organs and deep soft-tissue structures. Famous for its distinctive sweetheart heart-shaped distal paddle, this instrument distributes counter-traction evenly across a broad surface area, protecting fragile abdominal organs such as the liver, spleen, and intestines from focal compression injury.

Built for variable anatomical depths and surgical approaches, the Harrington Sweetheart Retractor is offered in five modular dimensions: 4 × 12.8 cm (\$129), 6.5 × 12.8 cm (\$134), 6.5 × 14.6 cm (\$139), 9 × 25.5 cm (\$144), and an extra-deep 11.5 × 25.5 cm (\$149) configuration. Surgeons can choose between an ergonomic curved shank that arches around wound margins or a rigid straight profile, both integrating seamlessly into modular table-mounted frames or manual handheld setups.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker manufactures the Harrington Sweetheart Retractor to withstand rigorous hospital autoclaving

cycles, providing general, urological, and gynecological surgical teams with long-lasting clinical durability and tactile control.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)