



RETRACTOR MAKER

PRODUCT BROCHURE OF

Dumbach Mandibular Channel Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Dumbach Mandibular Channel Retractor

SKU RM-DMCR-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Curved Left — 23 cm	RM-DMCR-1-1	\$49
Curved Right — 23 cm	RM-DMCR-1-2	\$49

PRODUCT DESCRIPTION

The Dumbach Mandibular Channel Retractor is a specialized maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic exposure of the inferior and posterior borders of the mandibular ramus during orthognathic and trauma procedures. Designed with a curved channel working end that hooks securely along the cortical bone border, this instrument shields facial vessels and the marginal mandibular nerve from rotating burrs and oscillating saws.

Built for deep intraoral clearance, this Dumbach Mandibular Channel Retractor configuration features an extended 23 cm total structural length (\$49), representing an advanced, lengthened design profile built to access deep mandibular fields comfortably while keeping the assistant's hands safely outside the operative corridor. Surgeons can customize the operational orientation at checkout between Curved Left and Curved Right to align with localized surgical approaches. Its ergonomic handle ensures steady, fatigue-free manual leverage.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Dumbach Mandibular Channel Retractor to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and plastic surgeons worldwide with superior durability, high corrosion resistance, and complete autoclave compatibility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



LinkedIn



Facebook



Instagram



X (Twitter)



YouTube



Pinterest