



RETRACTOR MAKER

PRODUCT BROCHURE OF

Fan Retractor With 3 Fingers

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Fan Retractor With 3 Fingers

SKU RM-FR3F-299-1

\$299 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
5 mm — Articulating	RM-FR3F-5A-1	\$299
5 mm — Non-Articulating	RM-FR3F-5N-1	\$299
10 mm — Articulating	RM-FR3F-10A-1	\$299
10 mm — Non-Articulating	RM-FR3F-10N-1	\$299

PRODUCT DESCRIPTION

The Fan Retractor With 3 Fingers is a specialized laparoscopic surgical retractor engineered specifically to provide dependable, focused soft-tissue displacement and visceral elevation during minimally invasive abdominal and pelvic procedures. Designed with three interlocking mechanical blades that close seamlessly into a low-profile cylinder, this instrument passes easily through standard trocars before fanning out into an anatomical triad that cradles tissues cleanly without crowding the optical field.

Built to clinical framework specifications, this Fan Retractor With 3 Fingers model (\$299) is offered in 5 mm or 10 mm shaft diameter options, with working lengths spanning from 330 mm to 375 mm (overall instrument length 438 mm to 487 mm). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument is selectable in Articulating or Non-Articulating shaft configurations, allowing surgeons to bend the working end around tight mesenteric angles. A dedicated luer-lock flush port is integrated into the handle assembly for pressurized lumen irrigation, while the smooth atraumatic finger borders prevent serosal abrasions or intestinal tears.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Fan Retractor With 3 Fingers to withstand repetitive hospital autoclave cycles, providing laparoscopic operating rooms, digestive surgery units,

and minimally invasive gynecological centers worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)