



RETRACTOR MAKER

PRODUCT BROCHURE OF

Obwegeser Ramus Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Obwegeser Ramus Retractor

SKU RM-ORR-22-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Obwegeser Ramus Retractor	22 cm overall length, anatomical curved ramus blade	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Obwegeser Ramus Retractor is an indispensable maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and osseous exposure along the ascending ramus of the mandible during orthognathic surgery. Designed with an anatomical hooking curve that seats securely over the posterior mandibular border, this instrument allows surgeons to perform bone cuts with reciprocating saws while shielding deep vascular structures and facial nerve branches from accidental trauma.

Built for uncompromised intraoral clearance, this Obwegeser Ramus Retractor configuration features an overall tool length of 22 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure lasting strength and corrosion resistance, its smoothly polished blade margins cradle the bone cortex cleanly, preventing mucosal tears or localized crush injury to the masseter and medial pterygoid muscles. The ergonomic handle provides operating assistants with steady manual leverage throughout extended osteotomies.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Obwegeser Ramus Retractor to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and facial reconstructive surgical suites worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)