



RETRACTOR MAKER

PRODUCT BROCHURE OF

Strandell Stille Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Strandell Stille Retractor

SKU RM-SSR-17-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Strandell Stille Retractor	17 cm overall length, 6 × 15 mm curved blade (Sharp or Blunt)	\$44
Protective Sleeve / Box	Individual transit packaging	\$44

PRODUCT DESCRIPTION

The Strandell Stille Retractor is a specialized handheld surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and wound margin exposure during plastic, hand, and general surgical interventions. Designed with a curved, tissue-conforming blade, this instrument allows surgical assistants to apply steady manual counter-traction across delicate dermal margins and fascial layers without causing localized pressure necrosis.

Built for high performance and lasting corrosion resistance, this dedicated Strandell Stille Retractor configuration features an overall tool length of 17 cm (\$44). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry layout incorporates a wide, flat blade with a curved, tissue-conforming configuration measuring 6 mm wide × 15 mm long. The available tissue interfaces include Sharp and Blunt design parameters, incorporating precise toothed profiles to deliver a secure, reliable hold on incision edges without slipping. Its ergonomic flat handle delivers fatigue-free manual grip stability.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Strandell Stille Retractor to withstand repetitive hospital autoclave cycles, providing plastic surgery suites, hand trauma centers, and general operating theaters worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)