



RETRACTOR MAKER

PRODUCT BROCHURE OF

Castaneda Rib Spreader

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Castaneda Rib Spreader

SKU RM-CRS-SERIES-1

\$309 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Option 1 — 10.5 mm × 32 mm — 65 mm maximum spread — Neonatal thoracotomy, coarctation repair, infant PDA ligation	RM-CRS-SERIES-1-1	\$309
Option 2 — 12.0 mm × 45 mm — 70 mm maximum spread — Infant congenital cardiac repair, vascular ring division	RM-CRS-SERIES-1-2	\$329
Option 3 — 15.0 mm × 60 mm — 120 mm maximum spread — Pediatric and adolescent thoracotomy, pulmonary shunt placement	RM-CRS-SERIES-1-3	\$349

PRODUCT DESCRIPTION

The Castaneda Rib Spreader is a specialized pediatric and neonatal thoracic surgical retractor engineered specifically to provide safe, controlled intercostal rib separation during congenital heart surgery and pediatric thoracotomies. Designed with delicate, finely calibrated slide tracks, this instrument allows pediatric cardiothoracic surgeons to distribute spreading forces with extreme micro-precision, avoiding rib fractures or costochondral dislocations in fragile infant thoracic cages.

Built to accommodate varying pediatric patient anatomies, this Castaneda Rib Spreader series is offered across three distinct size options: Option 1 (\$309) features 10.5 mm × 32 mm blades with a 65 mm maximum track spread for neonates; Option 2 (\$329) incorporates 12 mm × 45 mm blades with a 70 mm maximum track spread for infants; and Option 3 (\$349) features 15 mm × 60 mm blades with a 120 mm maximum track spread for pediatric and adolescent cases. Handcrafted from Premium Grade AISI 304/316 Stainless Steel

(Non-Magnetic), the smoothly curved lateral blades cradle the rib margins atraumatically, while the rotary screw drive locks firmly without back-slipping.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Castaneda Rib Spreader to withstand repetitive hospital sterilization cycles, providing pediatric cardiac operating rooms, children's hospitals, and congenital heart centers worldwide with exceptional structural durability, corrosion resistance, and complete autoclave dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)