



RETRACTOR MAKER

PRODUCT BROCHURE OF

Balfour USA Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Balfour USA Self Retaining Retractor

SKU RM-BUSRR-1

\$399 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
4 Sets (8 blades total)	Lateral Blades	65 mm, 85 mm, 105 mm, and 155 mm deep
1	Balfour Centre Blade	70 mm wide × 70 mm deep
1	Straight Centre Blade	70 mm wide × 100 mm deep
1	Balfour Frame	Central structural support frame with locking thumbscrew

PRODUCT DESCRIPTION

The Balfour USA Self Retaining Retractor is an internationally recognized abdominal surgical retractor platform engineered to provide reliable, hands-free soft-tissue exposure across exploratory laparotomies and pelvic operations. By utilizing a rigid rack-and-pinion sliding frame, this instrument locks bilateral wound edges open with mechanical counter-traction, freeing surgical assistants and maintaining stable exposure.

Designed for comprehensive procedural flexibility, this Balfour USA Self Retractor set features an interchangeable modular architecture. The system includes a solid Balfour frame, four pairs of lateral blades (65 mm, 85 mm, 105 mm, and 155 mm deep) to handle both slender and bariatric abdominal walls, a 70 × 70 mm Balfour center blade, and a 70 × 100 mm straight center blade for deep bladder and visceral retraction.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish, Retractor Maker supplies the Balfour USA Self Retaining Retractor non-sterile and fully autoclavable inside an individual protective box, providing surgical centers with rugged reliability and complete autoclave compatibility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)