



RETRACTOR MAKER

PRODUCT BROCHURE OF

Mitral Valve Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Mitral Valve Retractor

SKU RM-MVR-1

\$899 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Main Bracket	Main structural support component of the retractor system
1	Cooley Sternum Retractor	Sternum retraction component
1	Side Arm Bar	Provides additional structural support and positioning
3	Universal Clamp	Connects and secures compatible components within the system
1	Tricuspid Blade	Blade option for cardiac surgical exposure
1	Large Retractor Blade	Larger blade option for tissue retraction
3	Small Retractor Blade	Smaller blade options for tissue retraction
1	Sterilization Box	Organizes and protects the complete set during storage and reprocessing

PRODUCT DESCRIPTION

The Mitral Valve Retractor is a specialized cardiac surgery retractor system engineered to establish clear, stable exposure of the left atrium and intra-cardiac anatomy during valve reconstructions. By integrating directly onto a Cooley sternal spreader or sternotomy framework, this assembly elevates and isolates atrial tissue to give cardiac surgeons an unobstructed view of the valve annulus.

Configured for complete intraoperative positioning, the Mitral Valve Retractor system combines a heavy-duty Cooley sternum retractor, a primary structural support bracket, an adjustable side arm bar, and three universal mounting clamps. The blade set features a dedicated tricuspid blade, one large retractor blade, and three small retractor blades designed to gently retract atrial margins and interatrial septal tissue without causing tears.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-reflective satin finish, Retractor Maker supplies the complete Mitral Valve Retractor non-sterile and fully autoclavable in an organized, heavy-duty sterilization box, providing hospitals with consistent quality and dependable performance during cardiopulmonary bypass procedures.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)