



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hargis Mandibular Body Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hargis Mandibular Body Retractor

SKU RM-HMBR-17-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Hargis Mandibular Body Retractor	17 cm overall length, 13 mm hook tip — 1	\$44
Protective Sleeve / Box	Individual transit packaging — 1	\$44

PRODUCT DESCRIPTION

The Hargis Mandibular Body Retractor is an indispensable maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement along the mandibular corpus and lower border. Designed with an anatomical hook curve, this instrument allows surgeons to elevate periosteum and retract facial musculature smoothly, shielding the facial artery, anterior facial vein, and marginal mandibular nerve during bone division.

Built for precision operational control, this Hargis Mandibular Body Retractor configuration features an overall tool length of 17 cm (\$44). Its working end incorporates 13 mm tips featuring a versatile skin hook profile that seats securely along the cortical bone of the mandibular base, preventing tissue slippage without lacerating delicate periosteal tissues. The balanced ergonomic handle provides operating assistants with steady manual counter-traction, reducing hand fatigue during lengthy reconstructions.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Hargis Mandibular Body Retractor to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for oral, maxillofacial, and plastic surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)