



RETRACTOR MAKER

PRODUCT BROCHURE OF

Mathieu Finger Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Mathieu Finger Retractor

SKU RM-MFR-16-1

\$34 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Mathieu Finger Retractor	Dimensions / Details: 16 cm overall length, 8.0 mm width, 3 prongs each end (Sharp or Blunt) · Quantity Included: 1	\$34
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$34

PRODUCT DESCRIPTION

The Mathieu Finger Retractor is an essential double-ended surgical rake retractor engineered specifically to provide delicate, non-slip skin flap elevation and soft-tissue retraction during hand, plastic, and minor general surgical procedures. Designed with a balanced central handle connecting dual 3-prong working ends, this instrument allows surgical teams to transition effortlessly between differing tissue vectors without exchanging tools.

Built for optimal tactile control in confined operative fields, this Mathieu Finger Retractor configuration features an overall tool length of 16 cm (\$34) paired with an 8.0 mm working blade width. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), each working end incorporates 3 prongs customizable across selectable tip configurations: Blunt tips for gentle tissue manipulation minimizing local trauma, or Sharp tips for highly secure and firm tissue engagement in dense dermis. Its flat central handle ensures steady manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Mathieu Finger Retractor to withstand repetitive hospital sterilization cycles, providing outpatient surgical centers, plastic surgery suites, and hand trauma clinics worldwide with exceptional durability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)