



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sigmoid Notch Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sigmoid Notch Retractor

SKU RM-SNR-22-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Sigmoid Notch Retractor	22 cm overall length, 9.5 / 16 mm concave blade, fully bent tip	\$64
Protective Sleeve / Box	Individual transit packaging	\$64

PRODUCT DESCRIPTION

The Sigmoid Notch Retractor is a specialized oral and maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and osseous stabilization along the mandibular notch during ramus osteotomies. Designed to hook over the superior concavity between the coronoid process and the mandibular condyle, this instrument allows surgeons to perform vertical ramus or sagittal split procedures while completely shielding deep neurovascular structures from reciprocating saws.

Built to provide the perfect balance between clinical reach and manual control, this Sigmoid Notch Retractor configuration features an overall tool length of 22 cm (\$64). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure lasting strength and corrosion resistance, the working blade profile incorporates a 9.5 mm / 16 mm concave blade layout featuring a fully bent tip geometry that seats securely into the sigmoid notch cortex without slipping. The ergonomic flat handle provides operating assistants with steady manual leverage, preventing hand fatigue during prolonged midface and mandibular reconstructions.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Sigmoid Notch Retractor to withstand high-frequency autoclave sterilization protocols while offering a lifetime warranty against rust and defects, providing oral, maxillofacial, and facial reconstructive surgical teams worldwide with exceptional structural resilience, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)