



RETRACTOR MAKER

PRODUCT BROCHURE OF

EMORY Endo-Plastic Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



EMORY Endo-Plastic Retractor

SKU RM-EEPR-1

\$309 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
EMORY Endo-Plastic Retractor	10" (254 mm) length, 1" (25.4 mm) tip, 10 mm scope channel	\$309
Protective Sleeve / Box	Individual transit packaging	\$309

PRODUCT DESCRIPTION

The EMORY Endo-Plastic Retractor is a premier endoscopic aesthetic surgery instrument engineered specifically to provide expansive optical cavity elevation while securing a rigid optical corridor for endoscopic visualization. Designed for minimally invasive aesthetic dissections, this instrument allows surgeons to perform blunt subpectoral and subgaleal plane separation while keeping endoscopic camera optics clean and positioned directly over the dissection plane.

Built to rigorous surgical tolerances, this EMORY Endo-Plastic Retractor configuration features an extended 10-inch (254 mm) overall blade length (\$309) paired with a 1-inch (25.4 mm) wide blade tip. The instrument is engineered explicitly with an integrated scope guide channel optimized for 10 mm endoscopes, eliminating camera wobble and maintaining a stable, illuminated optical tunnel. Its ergonomic handle delivers refined tactile feedback, allowing surgeons to dissect tissue planes effortlessly without causing muscle tearing.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the EMORY Endo-Plastic Retractor non-sterile and fully autoclavable inside protective packaging, providing cosmetic, plastic, and reconstructive surgery centers with high durability, precision optics compatibility, and clinical longevity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)