



RETRACTOR MAKER

PRODUCT BROCHURE OF

Langenbeck Retractor-Solid Handle

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Langenbeck Retractor-Solid Handle

SKU RM-LRSH-203-1

\$49 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Style #1	Dimensions (Length × Width): 32 mm long × 10 mm wide · Primary Clinical Role: Focused wound margin and superficial soft-tissue retraction · Price (USD): \$49.00
1	Style #2	Dimensions (Length × Width): 38 mm long × 10 mm wide · Primary Clinical Role: Intermediate fascial retraction and tendon exposure · Price (USD): \$49.00
1	Style #3	Dimensions (Length × Width): 38 mm long × 12 mm wide · Primary Clinical Role: Standard adult appendectomy, hernia, and orthopedic joint access · Price (USD): \$49.00
1	Style #4	Dimensions (Length × Width): 38 mm long × 16 mm wide · Primary Clinical Role: Broad muscular displacement and deep soft-tissue holding · Price (USD): \$49.00
1	Protective Sleeve / Box	Dimensions (Length × Width): Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included

PRODUCT DESCRIPTION

The Langenbeck Retractor-Solid Handle is an essential handheld surgical retractor engineered specifically to provide dependable, non-slip soft-tissue and muscular retraction during general, orthopedic, and plastic surgery. Designed with an elongated right-angled blade terminating in a downward retaining lip, this instrument

allows surgical assistants to maintain constant counter-traction on cutaneous margins and muscle layers without slipping.

Built to clinical framework standards, this Langenbeck Retractor-Solid Handle features a uniform overall tool length of 20.3 cm (8 inches, \$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument is highly customizable through four distinct blade style classifications: Style #1 (32 mm long × 10 mm wide), Style #2 (38 mm long × 10 mm wide), Style #3 (38 mm long × 12 mm wide), and Style #4 (38 mm long × 16 mm wide). The solid, ergonomically balanced flat handle provides a comfortable, fatigue-free grip for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Langenbeck Retractor-Solid Handle to withstand repetitive hospital autoclave cycles, delivering high corrosion resistance, structural rigidity, and complete clinical durability for operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)