



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sims Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sims Vaginal Retractor

SKU RM-SVR-19-1

\$91 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
63.5 mm × 25 mm — Polished or Ebonized	RM-SVR-19-1-1	\$91
76.0 mm × 32 mm — Polished or Ebonized	RM-SVR-19-1-2	\$93
82.5 mm × 32 mm — Polished or Ebonized	RM-SVR-19-1-3	\$95
82.5 mm × 38 mm — Polished or Ebonized	RM-SVR-19-1-4	\$97
89.0 mm × 41 mm — Polished or Ebonized	RM-SVR-19-1-5	\$99

PRODUCT DESCRIPTION

The Sims Vaginal Retractor is a foundational gynecological and pelvic reconstructive surgical instrument engineered specifically to provide safe, non-traumatic displacement of the vaginal walls and pelvic floor. Designed with an asymmetric, double-ended layout consisting of two right-angled duckbill blades of different lengths facing the same direction, this instrument allows one end to be inserted deep into the vaginal vault while the opposing end serves as an ergonomic handle for manual traction.

Built to offer a lifetime warranty against rust and defects, this Sims Vaginal Retractor maintains an overall instrument length of 19 cm. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers two coating styles: a Standard Stainless Steel highly polished finish for routine procedures, or a non-reflective Ebonized Finish optimized for anti-glare laser surgical setups. The purchasing menu provides

five selectable blade dimensions: 63.5 mm × 25 mm (\$91), 76 mm × 32 mm (\$93), 82.5 mm × 32 mm (\$95), 82.5 mm × 38 mm (\$97), and 89 mm × 41 mm (\$99). The smooth, rounded blade margins protect mucosal folds from laceration.

Forged with exceptional precision, Retractor Maker manufactures the Sims Vaginal Retractor to withstand repetitive hospital autoclave cycles, providing obstetrical suites, gynecological theaters, and women's health clinics worldwide with unmatched durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)