



RETRACTOR MAKER

PRODUCT BROCHURE OF

Mollison Anterior Pillar Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Mollison Anterior Pillar Retractor

SKU RM-MAPR-22-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Mollison Anterior Pillar Retractor	22 cm overall length, dual 10 mm slender blades, angled handle	\$44
Protective Sleeve / Box	Individual transit packaging	\$44

PRODUCT DESCRIPTION

The Mollison Anterior Pillar Retractor is an essential otolaryngology surgical retractor engineered specifically to provide safe, non-traumatic displacement of the palatoglossal arch and anterior tonsillar pillar during tonsillectomy and oral procedures. Designed with an angled working configuration, this retractor allows ENT surgeons to displace mucosal margins cleanly, establishing an unencumbered visual corridor to the tonsillar bed for dissection and hemostasis.

Built to maximize precision and lower user hand fatigue, this Mollison Anterior Pillar Retractor configuration features an overall tool length of 22 cm (\$44). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates a blade width parameter of dual 10 mm slender blades contoured with rounded atraumatic margins to prevent mucosal tearing. The lightweight, angled center handle assembly provides surgical assistants with steady, fatigue-free manual leverage throughout extended operative cases.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Mollison Anterior Pillar Retractor to withstand repetitive hospital sterilization cycles, providing ENT operating suites, outpatient surgical clinics, and oral surgery centers worldwide with exceptional durability, corrosion resistance, and complete autoclave dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)