



RETRACTOR MAKER

PRODUCT BROCHURE OF

Army-Navy Retractor-Set of Two

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Army-Navy Retractor-Set of Two

SKU RM-ANR-SET-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Army-Navy Retractor	21 cm overall length, double-ended blades	\$49
Protective Sleeve / Box	Standard transit packaging	\$49

PRODUCT DESCRIPTION

The Army-Navy Retractor-Set of Two is a classic pair of double-ended surgical retractors engineered to provide steady, atraumatic tissue retraction during general, orthopedic, and plastic surgical procedures. By providing contrasting blade depths on either end of each instrument, this retractor pair allows surgical teams to transition effortlessly between shallow subcutaneous tissue holding and deeper muscular elevation.

Built for versatile clinical utility, this standard Army-Navy Retractor-Set of Two configuration features an overall length of 21 cm. Each instrument incorporates a central fenestrated handle for lightened weight and steady grip control, with one end presenting a short right-angled blade for superficial margins and the opposing end offering a longer right-angled blade with a slight lip to secure deep fascial margins without tissue slippage.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the Army-Navy Retractor-Set of Two non-sterile and fully autoclavable inside protective packaging, delivering high corrosion resistance and dependable durability for operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)