



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Shea Retractor With Fiber Optic Guide

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Shea Retractor With Fiber Optic Guide

**SKU** RM-SRFOG-22-1

**\$139 USD**

## COMPONENT LIST

| QUANTITY | COMPONENT                             | DIMENSIONS / DESCRIPTION                                  |
|----------|---------------------------------------|-----------------------------------------------------------|
| 1        | Shea Retractor With Fiber Optic Guide | 22 cm length, 14.5 mm wide blade, integrated light bundle |
| 1        | Universal Light Cable Adapter         | Connects to standard OR fiber-optic cables                |
| 1        | Protective Sleeve / Box               | Standard transit packaging                                |

## PRODUCT DESCRIPTION

The Shea Retractor With Fiber Optic Guide is a premier illuminated surgical retractor engineered specifically to combine non-traumatic soft-tissue displacement with direct, shadowless optical illumination in confined surgical corridors. Designed with an integrated fiber-optic channel along its shaft, this instrument directs intense, cool light straight to the working end, eliminating shadows in deep subfascial and facial planes without requiring secondary light wands.

Built to withstand frequent autoclave sterilization cycles while remaining highly resistant to corrosion, this Shea Retractor With Fiber Optic Guide configuration features an overall instrument length of 22 cm (\$139). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) paired with high-transmission optical glass fibers, the working blade width incorporates an expansive 14.5 mm wide blade layout contoured with smooth, rounded margins to protect delicate cutaneous flaps and nerves. The integrated lighting assembly features a universal connector compatible with standard OR light cables (Storz, Olympus, Wolf, ACMI).

Forged with an anti-glare satin matte finish on exposed metallic surfaces, Retractor Maker manufactures the Shea Retractor With Fiber Optic Guide to uphold rigorous clinical standards, providing plastic, reconstructive, and otolaryngology surgical centers worldwide with unmatched optical clarity, structural rigidity, and complete autoclave durability.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)