



RETRACTOR MAKER

PRODUCT BROCHURE OF

Levasseur Merrill Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Levasseur Merrill Retractor

SKU RM-LMR-25-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Levasseur Merrill Retractor	Dimensions / Details: 25 cm overall length, 12.5 / 27 mm osteotomy channel blade · Quantity Included: 1	\$69
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$69

PRODUCT DESCRIPTION

The Levasseur Merrill Retractor is an indispensable oral and maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and osseous exposure during intraoral vertical ramus osteotomies (IVRO). Designed with an anatomical hooking channel, this instrument allows surgeons to wrap securely around the posterior border of the ascending ramus, shielding the facial nerve branches, maxillary artery, and parotid tissue from reciprocating osteotomy saws.

Built for deep intraoral clearance, this Levasseur Merrill Retractor configuration features an overall tool length of 25 cm (\$69) that easily reaches deep mandibular fields while keeping assistant hands safely outside the operative field. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates specialized osteotomy parameter configurations featuring 12.5 mm / 27 mm dimension metrics optimized for vertical ramus osteotomies. Its ergonomic flat handle provides steady manual leverage throughout demanding bone cuts.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Levasseur Merrill Retractor to withstand repetitive hospital sterilization cycles, providing maxillofacial, orthognathic, and facial reconstructive surgeons worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)