



RETRACTOR MAKER

PRODUCT BROCHURE OF

Mollison Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Mollison Retractor

SKU RM-MLR-155-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Mollison Retractor	15.5 cm overall length, 4×4 sharp prongs, ratchet lock	\$99
Protective Sleeve / Box	Individual transit packaging	\$99

PRODUCT DESCRIPTION

The Mollison Retractor is an internationally recognized self-retaining surgical retractor engineered specifically to provide reliable, hands-free soft-tissue and periosteal retraction during mastoidectomy, thyroidectomy, and cranial procedures. Designed with curved spreading arms, this instrument holds incision margins open under constant mechanical tension, eliminating assistant hand fatigue and keeping the visual field unobstructed for high-speed drilling.

Built in a standard size layout engineered to provide a balanced distribution of mechanical pressure, this Mollison Retractor configuration features an overall tool length of 15.5 cm (\$99). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the prong layout configuration incorporates an asymmetrical 4 sharp prongs on the left and 4 sharp prongs on the right (4×4 array). The sharp prongs are meticulously crafted to grasp skin borders firmly and accurately, while the balanced 4×4 tension spreads pressure evenly to protect surrounding soft tissues from localized compression trauma. The base mechanical ratchet catch locks smoothly under counter-traction.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Mollison Retractor to withstand repetitive hospital autoclave cycles, providing ENT suites, neurosurgical operating theaters, and general surgical clinics worldwide with exceptional corrosion resistance, mechanical stability, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)