



RETRACTOR MAKER

PRODUCT BROCHURE OF

O'Sullivan-O'Conner Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



O'Sullivan-O'Conner Vaginal Retractor

SKU RM-OSVR-28-1

\$399 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Main Retractor Frame Chassis	Dimensions (Length × Width): 28 cm frame (Curved or Oval layout) · Quantity Included: 1 Frame · Price (USD): Included in Set (\$399.00)
1	Small Interchangeable Blade	Dimensions (Length × Width): 1.3 cm × 4.7 cm blunt blade · Quantity Included: 1 Blade · Price (USD): Included in Set
1	Medium Interchangeable Blade	Dimensions (Length × Width): 3.8 cm × 8.9 cm blunt blade · Quantity Included: 1 Blade · Price (USD): Included in Set
1	Large Interchangeable Blade	Dimensions (Length × Width): 4.4 cm × 8.9 cm blunt blade · Quantity Included: 1 Blade · Price (USD): Included in Set
1	Protective Storage / Transit Box	Dimensions (Length × Width): Standard Transit Packaging · Quantity Included: 1 Box · Price (USD): Included in Set

PRODUCT DESCRIPTION

The O'Sullivan-O'Conner Vaginal Retractor is a premier self-retaining surgical retractor system engineered specifically to provide expansive, hands-free soft-tissue displacement during complex gynecological, vaginal, and pelvic procedures. Designed with an adjustable perimeter chassis, this retractor maintains constant, multi-directional counter-traction against the pelvic wall and viscera, eliminating manual assistant fatigue and preserving a wide, static operative corridor.

Built to accommodate varying structural tissue depths, this O'Sullivan-O'Conner Vaginal Retractor set features a standardized overall frame dimension of 28 cm total length (\$399). The chassis frame geometry is available in either Curved or Oval configurations to adapt to diverse pelvic anatomies. The system includes three interchangeable blade parameters: a Small Blade (1.3 cm × 4.7 cm), a Medium Blade (3.8 cm × 8.9 cm), and a Large Blade (4.4 cm × 8.9 cm). All blades are built with blunt tips to guarantee safe, atraumatic soft tissue retraction without risk of mucosal perforation.

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker manufactures the O'Sullivan-O'Conner Vaginal Retractor to withstand repetitive hospital autoclave cycles, delivering exceptional mechanical stability, corrosion resistance, and complete clinical dependability for gynecological and urogynecological surgical centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)