



RETRACTOR MAKER

PRODUCT BROCHURE OF

Balfour Baby Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Balfour Baby Self Retaining Retractor

SKU RM-BBSRR-1

\$299 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1 Set (2 blades)	Lateral Blades	60 mm deep
1	Balfour Centre Blade	21 mm wide × 24 mm deep
1	Balfour Retractor Frame	125 × 90 mm

PRODUCT DESCRIPTION

The Balfour Baby Self Retaining Retractor is a specialized pediatric surgical retractor engineered to provide atraumatic, self-retaining exposure in neonatal, pediatric, and minor abdominal surgical fields. Built with scaled-down proportions, this instrument delivers steady, hands-free counter-traction without crowding delicate anatomical workspaces.

Designed specifically for smaller anatomical spaces, this complete Balfour Baby Self Retaining Retractor set features an agile 125 × 90 mm spreading frame, one pair of 60 mm deep lateral blades, and a compact 21 × 24 mm central blade for targeted pelvic and visceral displacement. The frame's fine-threaded thumb screw allows surgical teams to control spreading tension smoothly in millimeter increments, preventing tissue overstretching.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a glare-resistant satin matte finish, Retractor Maker supplies the Balfour Baby Self Retaining Retractor non-sterile and fully autoclavable inside an individual protective box, delivering exceptional clinical durability for pediatric surgery departments.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)