



RETRACTOR MAKER

PRODUCT BROCHURE OF

Petri Pterygoid Intra Oral Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Petri Pterygoid Intra Oral Retractor

SKU RM-PPIOR-23-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Petri Pterygoid Intra Oral Retractor	23 cm overall length, 25.5 mm wide pterygoid blade	\$79
Protective Sleeve / Box	Individual transit packaging	\$79

PRODUCT DESCRIPTION

The Petri Pterygoid Intra Oral Retractor is an essential oral, maxillofacial, and orthognathic surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and neurovascular shielding during Le Fort osteotomies. Designed to slide behind the maxillary tuberosity, this instrument shields the internal maxillary artery and pterygoid venous plexus from osteotomes during pterygomaxillary dysjunction.

Built to withstand heavy physical demands while maintaining handling ease, this Petri Pterygoid Intra Oral Retractor configuration features an overall tool length of 23 cm (\$79). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade width incorporates an expansive 25.5 mm wide blade layout optimized for Le Fort osteotome procedures. Its smoothly contoured blade edges wrap around the posterior maxillary wall cleanly, preventing soft-tissue lacerations while the ergonomic handle provides operating assistants with steady manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Petri Pterygoid Intra Oral Retractor to withstand repetitive hospital sterilization cycles, providing maxillofacial, orthognathic, and craniofacial surgery centers worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)