



RETRACTOR MAKER

PRODUCT BROCHURE OF

Lip And Cheek Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Lip And Cheek Retractor

SKU RM-LCR-26-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Lip And Cheek Retractor	Dimensions / Details: 26 cm overall length, curved atraumatic blade, flat handle · Quantity Included: 1	\$49
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$49

PRODUCT DESCRIPTION

The Lip And Cheek Retractor is an essential oral, maxillofacial, and dental surgical retractor engineered specifically to provide safe, non-traumatic displacement of the lips, cheeks, and oral mucosa during mandibular and orthognathic interventions. Designed with a curved blade that conforms smoothly against the buccal tissue envelope, this retractor shields the soft tissues from rotating rotary burrs and saws while maintaining an expansive visual corridor.

Built to optimize lower jaw accessibility in deep intraoral corridors, this dedicated Lip And Cheek Retractor configuration features an overall tool length of 26 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its extended reach allows operating assistants to maintain stable counter-traction from outside the surgical field, avoiding crowding around the dental arches. The smoothly polished blade margins distribute holding force evenly across mucosal layers, preventing localized tissue contusions or tears.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Lip And Cheek Retractor to withstand repetitive hospital sterilization cycles, providing dental, oral surgery, and facial reconstructive clinics worldwide with high corrosion resistance, unyielding structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)