



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sofield Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sofield Retractor

SKU RM-SFR-SERIES-1

\$46 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Size #1 — 18.0, 19.5, 20.3, 21.0, 22.3, 22.8 cm	RM-SFR-SERIES-1-1	\$46
Size #2 — 18.0, 19.5, 20.3, 21.0, 22.3, 22.8 cm	RM-SFR-SERIES-1-2	\$46
Size #3 — 18.0, 19.5, 20.3, 21.0, 22.3, 22.8 cm	RM-SFR-SERIES-1-3	\$46
Size #4 — 18.0, 19.5, 20.3, 21.0, 22.3, 22.8 cm	RM-SFR-SERIES-1-4	\$46
Size #5 — 18.0, 19.5, 20.3, 21.0, 22.3, 22.8 cm	RM-SFR-SERIES-1-5	\$46
Size #6 — 18.0, 19.5, 20.3, 21.0, 22.3, 22.8 cm	RM-SFR-SERIES-1-6	\$46

PRODUCT DESCRIPTION

The Sofield Retractor is an essential orthopedic surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue and periosteal displacement along long bone diaphyses during osteotomies and fracture repairs. Designed with an elongated curved blade that wraps smoothly around the bone cortex, this instrument allows surgeons to hold muscular bellies and periosteal flaps away from bone cuts without slipping or causing focal neurovascular compression.

Built to withstand frequent autoclave sterilization cycles while remaining highly resistant to corrosion, this Sofield Retractor line is highly customizable across six distinct numbered model sizes (\$46). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), all models (Size #1 through Size #6) are selectable in six overall lengths: 18.0 cm, 19.5 cm, 20.3 cm, 21.0 cm, 22.3 cm, and 22.8 cm. The active blade dimensions range from 13 mm × 19 mm in Size #1 up to an expansive 19 mm × 114 mm in Size #6. Each variant features an integrated flat surface comfort grip handle that delivers superior manual leverage and reduces hand fatigue during lengthy bone reconstructions.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Sofield Retractor to uphold strict clinical standards, providing pediatric orthopedics, trauma centers, and long-bone reconstructive suites worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)