



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Anterior Femoral Neck Anteromedial Rim Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Anterior Femoral Neck Anteromedial Rim Retractor

**SKU** RM-AFNAR-1

**\$49 USD**

## AVAILABLE SIZES & PRICING

| SIZE / OPTION                                    | SKU / DETAILS                                         | PRICE (USD) |
|--------------------------------------------------|-------------------------------------------------------|-------------|
| Anterior Femoral Neck Anteromedial Rim Retractor | 15" (37.46 cm) length, curved blade, hexagonal handle | \$49        |
| Protective Sleeve / Box                          | Individual transit packaging                          | \$49        |

## PRODUCT DESCRIPTION

The Anterior Femoral Neck Anteromedial Rim Retractor is a heavy-duty orthopedic surgical instrument engineered specifically to expose the anterior femoral neck and anteromedial acetabular rim during total hip replacement procedures. Designed with an anatomical curvature that conforms around the proximal femur, this instrument allows surgical assistants to maintain steady leverage while retracting soft tissue away from the femoral head and neck.

Engineered for deep-cavity pelvic and hip procedures, this Anterior Femoral Neck Anteromedial Rim Retractor configuration features an extended 15-inch (37.46 cm) blade scale that keeps the assistant's hands safely outside the surgeon's operative field. The curved and angled distal blade is contoured to seat securely against the anteromedial bony ridge without damaging adjacent neurovascular bundles. Its specialized hexagonal handle ensures non-slip rotational control, optimal grip stability, and reduced hand fatigue during prolonged retraction.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish, Retractor Maker manufactures the Anterior Femoral Neck Anteromedial Rim Retractor to endure heavy orthopaedic torque and repeated steam sterilization cycles, providing joint reconstruction teams with long-lasting clinical dependability.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)