



RETRACTOR MAKER

PRODUCT BROCHURE OF

Direct Access Anterior Broaching Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Direct Access Anterior Broaching Retractor

SKU RM-DAABR-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Left Orientation — 9" (22.86 cm), 57 mm blade	RM-DAABR-1-1	\$99
Right Orientation — 9" (22.86 cm), 57 mm blade	RM-DAABR-1-2	\$99

PRODUCT DESCRIPTION

The Direct Access Anterior Broaching Retractor is a specialized orthopedic surgical instrument engineered specifically to provide clear, stable soft-tissue and muscular clearance during the femoral preparation phase of direct anterior approach (DAA) total hip arthroplasty. Designed to elevate the proximal femur into the wound opening, this retractor allows surgeons to broach the femoral canal without encountering tissue impingement or damaging the tensor fasciae latae.

Built for precision joint leverage, this Direct Access Anterior Broaching Retractor configuration features a 9-inch (22.86 cm) total structural reach (\$99) paired with an expansive 57 mm wide blade profile. Its working end is formatted with a specialized curved pointed tip engineered to leverage bone contours smoothly, hooking securely behind the greater trochanter or femoral neck. Surgeons can choose between dedicated Left and Right operational orientations to align perfectly with the surgical approach.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Direct Access Anterior Broaching Retractor to withstand immense orthopedic torque and repeated hospital autoclave cycles, providing hip reconstruction surgeons worldwide with exceptional durability and clinical stability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)