



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sandhu Nerve Root Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sandhu Nerve Root Retractor

SKU RM-SNRR-243-1

\$54 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Sandhu Nerve Root Retractor	24.3 cm overall length, 3.2 mm blunt blade, 90° malleable shaft	\$54
Protective Sleeve / Box	Individual transit packaging	\$54

PRODUCT DESCRIPTION

The Sandhu Nerve Root Retractor is a specialized spinal microsurgery instrument engineered specifically to provide safe, non-traumatic retraction and protection of delicate nerve roots and the thecal sac during posterior lumbar discectomy and decompression. Designed with a slender malleable shaft, this retractor allows spine surgeons to bend and adjust the working angle intraoperatively to navigate challenging interlaminar corridors without obstructing the microscopic line of sight.

Built for deep operative reaching clearance and micro-precision, this Sandhu Nerve Root Retractor configuration features an overall tool length of 24.3 cm (\$54). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to provide exceptional durability and corrosion resistance, its working blade profile incorporates a 3.2 mm wide blunt blade layout engineered to ensure gentle tissue handling and minimize local trauma. The shaft geometry features a 90° angled malleable shaft configuration that holds its formed shape under tissue tension without springback. Its ergonomic handle delivers refined tactile feedback.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Sandhu Nerve Root Retractor to withstand repetitive hospital autoclave cycles, providing neurosurgical suites and orthopedic spine centers worldwide with unyielding structural reliability, high corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)