



RETRACTOR MAKER

PRODUCT BROCHURE OF

Lahey Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Lahey Retractor

SKU RM-LR-HH-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
18 cm (approx. 7")	RM-LR-HH-1-1	\$39
20.3 cm (approx. 8")	RM-LR-HH-1-2	\$49

PRODUCT DESCRIPTION

The Lahey Retractor is an internationally recognized handheld surgical retractor engineered to provide atraumatic soft-tissue exposure in delicate neck and endocrine procedures. Designed with a slender right-angled working tip, this instrument allows surgeons to gently elevate strap muscles, fascia, and cutaneous margins without obscuring visual access to underlying neurovascular structures.

Optimized for precision cervical dissections, the Lahey Retractor is available in two distinct overall working lengths: 18 cm (\$39) and 20.3 cm (\$49). The retractor features selectable blade configurations including a narrow 5 mm blade and a 15 mm × 3 mm right-angle blade, providing targeted holding capability with an ergonomic flat handle that ensures comfortable, fatigue-free control for operating assistants.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish, Retractor Maker manufactures the Lahey Retractor to strict clinical standards, delivering high corrosion resistance, unyielding rigidity, and full autoclave compatibility for general surgery and otolaryngology departments.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)