



RETRACTOR MAKER

PRODUCT BROCHURE OF

McCabe Antral Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



McCabe Antral Retractor

SKU RM-MAR-658-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard — 6 5/8 in	RM-MAR-658-1-1	\$99

PRODUCT DESCRIPTION

The McCabe Antral Retractor is an indispensable otolaryngology surgical retractor engineered specifically to provide safe, hands-free soft-tissue and cheek retraction during maxillary sinus surgery and Caldwell-Luc interventions. Designed with specialized curved shanks, this instrument elevates labial and cheek tissues cleanly away from the canine fossa, establishing an unobstructed operative corridor to the maxillary antrum.

Built for absolute structural control, this McCabe Antral Retractor configuration features an overall tool length of 6 5/8 inches (16.83 cm, \$99). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates curved shanks integrated with a 2-blade centralized retraction layout that holds mucosal margins firmly under continuous mechanical tension. The handle assembly features a finger ring grip paired with an integrated multi-position ratchet gear lock that eliminates assistant hand fatigue throughout extensive antral explorations.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the McCabe Antral Retractor to withstand repetitive hospital sterilization cycles, providing ENT suites, maxillofacial clinics, and hospital operating theaters worldwide with exceptional corrosion resistance, mechanical stability, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)