



RETRACTOR MAKER

PRODUCT BROCHURE OF

Minnesota Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Minnesota Retractor

SKU RM-MNR-55-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard Pattern (Without Suction)	RM-MNR-55-1-1	\$39
With Suction Tube	RM-MNR-55-1-2	\$49

PRODUCT DESCRIPTION

The Minnesota Retractor is an internationally acclaimed oral, dental, and maxillofacial surgical retractor engineered specifically to provide dependable, non-traumatic displacement of the cheeks, lips, and tongue during intraoral procedures. Designed with an offset double-angled shank, this instrument allows surgical assistants to shield delicate buccal tissues cleanly away from rotating dental burs and oscillating saws while maintaining an expansive visual corridor.

Built to clinical framework standards, this Minnesota Retractor maintains a uniform longitudinal profile of 5.5 inches (13.97 cm) paired with a 20.5 mm wide working blade. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers two distinct structural typologies selectable from the option menus: a Standard Pattern (Without Suction) featuring a solid-surface baseline model (\$39), or a configuration With Suction Tube (\$49) integrating a continuous fluid-evacuation channel directly onto the retractor blade. The ergonomic central handle ensures fatigue-free manual leverage and slip-resistant handling.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Minnesota Retractor to withstand repetitive hospital sterilization cycles, providing oral surgery suites, maxillofacial clinics, and dental surgical facilities worldwide with exceptional durability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)