



RETRACTOR MAKER

PRODUCT BROCHURE OF

Redundant Wall Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Redundant Wall Retractor

SKU RM-RWR-18-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Medium — 18 cm, 5.7×2.2 cm blade	RM-RWR-18-1-1	\$89
Large — 18 cm, 7.0×2.2 cm blade	RM-RWR-18-1-2	\$99
X-Large — 18 cm, 10.8×3.0 cm blade	RM-RWR-18-1-3	\$109

PRODUCT DESCRIPTION

The Redundant Wall Retractor is an essential gynecological and urogynecological surgical retractor engineered specifically to prevent redundant lateral vaginal tissue folds from prolapsing into the operative field. Designed with an elongated curved blade, this instrument allows surgeons to support and retract the lateral vaginal walls cleanly, establishing an unhindered visual corridor to the cervix and fornices during complex vaginal procedures.

Built in a dedicated configuration profile to promote a hygienic working surface, this Redundant Wall Retractor configuration features an overall tool length of 18 cm (\$89). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade width measures 2 cm wide, with a working blade length range scaled from 9 cm to 12.5 cm long to match variable vaginal capacities. Its smooth, highly polished finish minimizes mucosal friction, while the ergonomic flat handle provides operating assistants with steady manual leverage throughout prolonged pelvic floor reconstructions.

Forged with an anti-glare satin matte finish on handling surfaces, Retractor Maker produces the Redundant Wall Retractor to withstand repetitive hospital autoclave cycles, providing obstetrical, gynecological, and urogynecological surgical centers worldwide with exceptional durability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)