



RETRACTOR MAKER

PRODUCT BROCHURE OF

Surgical Lead Hand Retractor Plate

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Surgical Lead Hand Retractor Plate

SKU RM-SLHRP-199-1

\$199 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Surgical Lead Hand Retractor Plate	Adult hand-shaped malleable lead retractor plate	\$199
Protective Sleeve / Box	Standard transit packaging	\$199

PRODUCT DESCRIPTION

The Surgical Lead Hand Retractor Plate is an indispensable orthopedic and hand surgery positioning device engineered specifically to provide reliable, hands-free immobilization of the hand and fingers during intricate reconstructive and microvascular interventions. Designed with individual bendable finger tabs, this plate allows surgeons to contour the malleable fingers around the patient's digits, securing the hand in full extension, abduction, or flexion without requiring assistant holding.

Built for repeated clinical reshaping and long-life durability, this Surgical Lead Hand Retractor Plate configuration (\$199) provides a heavy, stable operating platform that rests flat on the hand table. Handcrafted from premium medical-grade malleable lead with a non-porous protective surface finish, the pliable finger extensions can be bent and unbent thousands of times without metal fatigue or cracking. The smooth, rounded edges prevent focal pressure necrosis on delicate palmar skin and dorsal digital vessels.

Manufactured to withstand repetitive hospital sterilization cycles, Retractor Maker produces the Surgical Lead Hand Retractor Plate to uphold strict clinical standards, providing hand surgery clinics, burn reconstructive units, and orthopedic trauma operating rooms worldwide with exceptional positioning stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)