



RETRACTOR MAKER

PRODUCT BROCHURE OF

Rultract Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Rultract Retractor System

SKU RM-RRS-1

\$1,499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Spline Tip, Double Blend Post/Clamp	Provides the primary mounting interface and clamp adjustment
1	Rake Assembly with Snap Clip	Engages and holds surgical tissue with a snap-clip attachment
1	Spline Top Post	Connects the upper frame section and stabilizes the post assembly

PRODUCT DESCRIPTION

The Rultract Retractor System is a specialized table-mounted surgical retractor engineered to provide consistent, hands-free upward tissue elevation and unmatched anatomical visibility during demanding operative procedures. By anchoring securely to standard operating table rails, this modular platform eliminates surgical assistant fatigue while maintaining smooth, non-slip vertical and lateral traction throughout prolonged interventions.

Specifically configured for sternal lifting, internal mammary artery (IMA) dissection, and upper abdominal exposures, this Rultract Retractor System includes a heavy-duty spline tip double blend post with integrated clamp, a stabilizing spline top post, and an articulating rake assembly featuring a secure snap clip. Each component is machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) treated with an anti-glare satin finish, preventing harsh reflections under surgical overhead lighting and safeguarding surrounding anatomical structures from excessive focal tension.

Supplied non-sterile and fully autoclavable, Retractor Maker manufactures the Rultract Retractor System to exact clinical tolerances, delivering high-performance retraction instruments directly to hospitals, cardiac centers, and surgical distributors worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)