



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hohmann Retractor Square End

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hohmann Retractor Square End

SKU RM-HRSE-1

\$74 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8 ½" (21.59 cm) — Curved blade with square distal end	RM-HRSE-1-1	\$74
9 ½" (24.13 cm) — Curved blade with square distal end	RM-HRSE-1-2	\$79

PRODUCT DESCRIPTION

The Hohmann Retractor Square End is a specialized orthopedic surgical instrument engineered specifically to provide broad, stable bone leverage and soft-tissue displacement during joint arthroplasty and long bone fracture repair. Designed with a distinct squared distal end, this retractor distributes counter-traction forces across a wide cortical contact surface, preventing localized bone indentation or slippage on flat osseous planes.

Built for reliable joint space extraction, this Hohmann Retractor Square End model is available in two distinct overall tool lengths: an 8 ½-inch (21.59 cm) configuration (\$74) for standard exposures and an extended 9 ½-inch (24.13 cm) configuration (\$79) for deeper operative fields. Its square-ended blade seats firmly against the femoral or tibial cortex, allowing surgical assistants to elevate heavy muscular planes with minimal effort while maintaining a clear, static visual corridor.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Hohmann Retractor Square End to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing orthopedic surgery centers worldwide with unyielding clinical durability, high corrosion resistance, and complete dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)