



RETRACTOR MAKER

PRODUCT BROCHURE OF

Levasseur Merrill Retractor Fiber Optic

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Levasseur Merrill Retractor Fiber Optic

SKU RM-LMR-FO-1

\$139 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Levasseur Merrill Retractor Fiber Optic	Dimensions / Details: 25 cm length, 12.5 / 27 mm blade, fiber-optic light bundle · Quantity Included: 1
1	Universal Light Cable Adapter	Dimensions / Details: Connects to standard OR light cables · Quantity Included: 1
1	Protective Sleeve / Box	Dimensions / Details: Standard transit packaging · Quantity Included: 1

PRODUCT DESCRIPTION

The Levasseur Merrill Retractor Fiber Optic is a premier illuminated maxillofacial surgical retractor engineered specifically to provide safe soft-tissue protection combined with direct, shadowless optical illumination during vertical ramus osteotomies. By channeling high-intensity light along the inner contour of its channel blade, this instrument eliminates shadowing in deep intraoral tunnels, allowing surgeons to execute vertical osteotomy cuts with absolute anatomical confidence.

Built for uncompromised mandibular bone interventions, this Levasseur Merrill Retractor Fiber Optic model features an overall instrument length of 25 cm (\$139). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates specialized osteotomy parameter configurations featuring 12.5 mm / 27 mm dimension metrics explicitly optimized for vertical osteotomies. Its integrated fiber-optic bundle includes a universal connector compatible with standard OR light cables (Storz, Olympus, Wolf, ACMI), ensuring effortless connectivity to light sources.

Crafted with an anti-glare satin matte finish on all metallic surfaces, Retractor Maker produces the Levasseur Merrill Retractor Fiber Optic to withstand repetitive hospital sterilization cycles, providing maxillofacial and orthognathic surgical centers worldwide with exceptional optical clarity, unyielding structural rigidity, and

complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)