



RETRACTOR MAKER

PRODUCT BROCHURE OF

Medial Ramus Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Medial Ramus Retractor

SKU RM-MRR-195-1

\$99 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Left Model Variant	Overall Length: 19.5 cm · Primary Clinical Role: Medial soft-tissue protection in left-sided sagittal split osteotomies · Price (USD): \$99.00
1	Right Model Variant	Overall Length: 19.5 cm · Primary Clinical Role: Medial soft-tissue protection in right-sided sagittal split osteotomies · Price (USD): \$99.00
1	Protective Sleeve / Box	Overall Length: Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included

PRODUCT DESCRIPTION

The Medial Ramus Retractor is an indispensable orthognathic surgical instrument engineered specifically to provide safe, non-traumatic soft-tissue displacement along the medial aspect of the ascending ramus during sagittal split procedures. Designed to slide gently between the medial ramus cortex and the medial pterygoid muscle, this retractor shields the inferior alveolar neurovascular bundle above the lingula during horizontal bone cuts.

Built to accommodate bilateral surgical approaches, this Medial Ramus Retractor model features an overall tool length of 19.5 cm (\$99). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers distinct anatomical variations selectable from the configuration menu: a Left model variant and a Right model variant. The anatomically curved blade tip seats snugly around the posterior and medial border of the ramus, preventing slippage while dropping assistant hands safely away from the surgeon's line of sight.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Medial Ramus Retractor to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and facial reconstructive surgeons worldwide with exceptional durability, corrosion resistance, and complete autoclave dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)