



RETRACTOR MAKER

PRODUCT BROCHURE OF

Contour Self-Retaining Scalp Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Contour Self-Retaining Scalp Retractor

SKU RM-CSRSR-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
3×4 Blunt Prongs	RM-CSRSR-1-1	\$99
3×4 Sharp Prongs	RM-CSRSR-1-2	\$99

PRODUCT DESCRIPTION

The Contour Self-Retaining Scalp Retractor is an essential neurosurgical instrument engineered specifically to provide steady, hands-free exposure of the cranium and underlying dural structures during craniotomies and skull base procedures. Designed with a straight body layout, this retractor supports precise structural alignment across flat incisions along the scalp, maintaining broad operative access without slippage.

Built for complete procedural adaptability, this Contour Self-Retaining Scalp Retractor model features an overall length of 14 cm (\$99) that provides optimal working space while preventing operative crowding. Surgeons can configure the working ends with either 3×4 Blunt Prongs—ideal for providing gentle yet effective tissue retraction near delicate superficial vessels—or 3×4 Sharp Prongs optimized to deliver maximum stability and strong leverage on dense scalp tissue structures. Its finger-ring ratchet catch locks firmly under tension, eliminating manual assistant strain.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Contour Self-Retaining Scalp Retractor to withstand repetitive hospital sterilization cycles, providing neurosurgical and spinal operating theaters with exceptional durability, smooth ratchet action, and clinical stability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)