



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Custom Bookwalter Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Custom Bookwalter Retractor

**SKU** RM-CBR-1

**\$2,999 USD**

## COMPONENT LIST

| QUANTITY     | COMPONENT                                    | DIMENSIONS / DESCRIPTION                                                                                                                                            |
|--------------|----------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Configurable | Mounting Hardware (available options)        | Standard Field Posts (12", 18"), Top Mount StrongArm PLUS Holders, Heavy-Duty Rail Clamps, Post Couplers, Horizontal Bars, Horizontal Flex Bars, Extension Bars     |
| Configurable | Ring Frames (available options)              | Solid Oval Rings, Solid Round Rings, Segmented Oval Rings, Segmented Round Rings, Half-Circle Segments, Straight Ring Segments (4", 6", 8"), Notched Perineal Rings |
| Configurable | Ratchet Mechanisms (available options)       | Standard Non-Tilt Ratchets, Multi-Angle Tilt Ratchets, Twist Ratchets, Ratchet Extenders, Low-Profile Lock Mechanisms                                               |
| Configurable | Rigid & Malleable Blades (available options) | Kelly, Balfour, Harrington, Deaver, Malleable, Richardson, Meyerding, Gelpi, Vaginal Lateral/Posterior, Renal Wylie and Rectal Blades, Deep Rake Claws              |
| Configurable | Storage & Packaging (available options)      | Perforated Sterilization Boxes, Customized Tray Inserts, Protective Transport Containers, OEM Private Label Packaging                                               |

## PRODUCT DESCRIPTION

The Custom Bookwalter Retractor is a fully configurable table-mounted self-retaining retractor system manufactured according to the customer's surgical specialty, preferred components and procedural requirements. Rather than a fixed set, it is specified together with you: you choose the mounting hardware, ring frame style, ratchet mechanisms, rigid and malleable blades, and the storage or packaging format.

Important: the component categories listed under Component List are the options available for selection. They are not all included at the \$2,999 reference price. The exact contents of your set — and the final price — are confirmed in writing during quoting, once the configuration is agreed.

Available mounting hardware includes standard 12" and 18" field posts, top-mount StrongArm PLUS holders, heavy-duty rail clamps, post couplers, horizontal bars, horizontal flex bars and extension bars. Ring frames can be solid or segmented, oval or round, and include half-circle segments, 4", 6" and 8" straight segments and notched perineal rings. Ratchet options cover standard non-tilt, multi-angle tilt and twist ratchets, ratchet extenders and low-profile lock mechanisms. Blade options include Kelly, Balfour, Harrington, Deaver, malleable, Richardson, Meyering, Gelpi, vaginal lateral and posterior, renal Wylie and rectal blades, plus deep rake claws.

All components are manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic). Sets are supplied non-sterile, are reusable and fully autoclavable, and must be cleaned and sterilized according to your facility's validated procedures before use. Storage options include perforated sterilization boxes, customized tray inserts, protective transport containers and OEM private-label packaging.

Retractor Maker manufactures custom Bookwalter-style sets directly at its facility near Sardaar Saleem Medident Hospital, Roras Road, Sialkot 51310, Pakistan.

Retractor Maker manufactures and supplies the Custom Bookwalter Retractor for hospitals, distributors and medical suppliers. The Custom Bookwalter Retractor is manufactured, inspected and supplied directly by Retractor Maker.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractor maker.com](mailto:info@retractor maker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractor maker.com](http://retractor maker.com)**  
**[retractor manufacturer.com](http://retractor manufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



LinkedIn



Facebook



Instagram



X (Twitter)



YouTube



Pinterest