



RETRACTOR MAKER

PRODUCT BROCHURE OF

Rochard Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Rochard Retractor

SKU RM-RR-1

\$1,699 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
2	Rail Clamp	Operating table rail attachment
2	Table Side Bar	26" (Mounting Frame)
1	Table Cross Bar	21" (Mounting Frame)
1	Rochard Fixation Device	Frame fixation mechanism
1	Rochard Abdominal Retractor	50 × 90 mm
1	Rochard Abdominal Retractor	50 × 105 mm
1	Rochard Abdominal Retractor	50 × 120 mm

PRODUCT DESCRIPTION

The Rochard Retractor is a classic, table-mounted self-retaining retraction frame engineered to deliver stable upward elevation and broad lateral exposure during upper abdominal procedures. Anchoring firmly to both sides of the surgical table, its bilateral bridge construction eliminates the continuous physical strain of manual retraction while maintaining rigid exposure of the subcostal margin and deep visceral structures.

Built for heavy abdominal wall retraction, the Rochard Retractor system consists of two robust rail clamps, dual 26-inch vertical/side bars, a 21-inch crossbar, and a dedicated Rochard fixation device. The mechanical screw-tensioning device allows surgeons to exert powerful, micro-adjustable upward leverage on the three included broad abdominal blades (50 × 90 mm, 50 × 105 mm, and 50 × 120 mm) to lift the costal arch without risk of slipping.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin finish, Retractor Maker manufactures the Rochard Retractor to exacting hospital specifications, providing surgical teams worldwide with dependable durability and complete autoclave compatibility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)