



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hibbs Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hibbs Retractor

SKU RM-HBR-24-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
51 mm long × 19 mm wide — 24 cm	RM-HBR-24-1-1	\$39
76 mm long × 25 mm wide — 24 cm	RM-HBR-24-1-2	\$41
82 mm long × 25 mm wide — 24 cm	RM-HBR-24-1-3	\$43
101 mm long × 25 mm wide — 24 cm	RM-HBR-24-1-4	\$45

PRODUCT DESCRIPTION

The Hibbs Retractor is a foundational orthopedic and spinal surgical retractor engineered specifically to provide reliable, deep soft-tissue displacement and paraspinous muscle exposure during spine and joint reconstructions. Designed with an elongated right-angled blade terminating in a specialized toothed retaining lip, this instrument anchors securely against dense bone cortices, preventing tissue slippage during vigorous bone reaming and decortication.

Built for heavy surgical counter-traction, this Hibbs Retractor series maintains a uniform instrument length of 24 cm. The retractor features an array of selectable working blade lengths and widths available individually or as a multi-tool package: 51 mm long × 19 mm wide (\$39), 76 mm long × 25 mm wide (\$41), 82 mm long × 25 mm wide (\$43), and an extra-deep 101 mm long × 25 mm wide (\$45) profile. Its teardrop-aperture ring handle provides surgical assistants with superior manual pulling leverage, reducing hand fatigue throughout lengthy procedures.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Hibbs Retractor to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing spinal and orthopedic trauma centers worldwide with unyielding durability,

corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)