



RETRACTOR MAKER

PRODUCT BROCHURE OF

Desmarres EyeLid Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Desmarres EyeLid Retractor

SKU RM-DER-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard Configuration — 5 ½" (14 cm) length, stainless steel saddle blade — Routine eyelid retraction and globe protection	RM-DER-1-1	\$49
Insulated Configuration — 5 ½" (14 cm) length, dielectric polymer coating — Electrosurgical safety during cautery and radiofrequency cutting	RM-DER-1-2	\$79

PRODUCT DESCRIPTION

The Desmarres EyeLid Retractor is an internationally acclaimed ophthalmic and plastic surgery retractor engineered specifically to provide safe, non-traumatic retraction of the upper or lower eyelid during aesthetic and reconstructive eyelid procedures. Designed with an ultra-smooth, saddle-shaped blade, this instrument cradles the tarsal plate and conjunctiva comfortably, shielding the delicate ocular globe from surgical trauma.

Built for refined periocular handling, this Desmarres EyeLid Retractor configuration features an overall length of 5 ½ inches (14 cm) that gives the surgeon exceptional fingertip mobility. The instrument is available in two specialized versions: a Standard Stainless Steel model (\$49) for routine dissection and an Insulated Configuration (\$79) coated with a durable dielectric polymer to safeguard patient tissues from accidental thermal transfer during monopolar or bipolar electrosurgery. Its balanced handle reduces hand fatigue during intricate orbital work.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker produces the Desmarres EyeLid Retractor to strict clinical standards, delivering high corrosion resistance, smooth rounded margins, and complete autoclave durability for cosmetic, oculoplastic,

and ophthalmic surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)