



RETRACTOR MAKER

PRODUCT BROCHURE OF

MIS Nerve Root Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



MIS Nerve Root Retractor

SKU RM-MIS-NRR-17-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
4.7 mm Wide Blade	RM-MIS-NRR-17-1-1	\$84
7.3 mm Wide Blade	RM-MIS-NRR-17-1-2	\$89

PRODUCT DESCRIPTION

The MIS Nerve Root Retractor is a specialized spinal microsurgical instrument engineered specifically to provide safe, non-traumatic nerve root and thecal sac displacement through minimally invasive tubular portals. Designed with an ultra-slender shaft and an angled distal tip, this retractor allows spine surgeons to manipulate delicate neural elements within narrow working corridors without obstructing microscopic illumination or surgical instruments.

Built for restricted spinal pathways, this MIS Nerve Root Retractor retains a uniform longitudinal profile of 17 cm that fits comfortably inside tubular dilator systems. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers two selectable blade width profiles: a 4.7 mm wide blade layout (\$84) for narrow micro-corridors and exiting nerve root isolation, and an expanded 7.3 mm wide blade layout (\$89) for broad traversing nerve root and thecal sac shielding. The smooth, rounded blade margins prevent localized traction neuropraxia or dural tears.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the MIS Nerve Root Retractor to withstand repetitive hospital sterilization cycles, providing neurosurgical and orthopedic spine centers worldwide with unyielding structural rigidity, high corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)