



RETRACTOR MAKER

PRODUCT BROCHURE OF

Miller Rake Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Miller Rake Retractor

SKU RM-MRR-514-1

\$39 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Blunt Tips Profile	Dimensions / Details: 5 ¼" length, 15 mm width, 9 blunt prongs · Primary Clinical Role: Safe cutaneous flap elevation with zero puncture risk · Price (USD): \$39.00
1	Sharp Tips Profile	Dimensions / Details: 5 ¼" length, 15 mm width, 9 sharp prongs · Primary Clinical Role: Secure non-slip purchase on dense dermis and fibrous fascia · Price (USD): \$39.00
1	Protective Sleeve / Box	Dimensions / Details: Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included

PRODUCT DESCRIPTION

The Miller Rake Retractor is a specialized handheld dermatologic and plastic surgery retractor engineered specifically to provide broad, non-slip skin flap elevation and subcutaneous tissue exposure. Designed with an expansive 9-prong rake head, this instrument allows surgeons to apply uniform manual counter-traction across cutaneous and subcutaneous planes without causing focal dermal tearing or margin ischemia.

Built in a standard size configuration to optimize manual dexterity, this Miller Rake Retractor configuration features an overall instrument length of 5 ¼ inches (13.3 cm) paired with a 15.0 mm working end width (\$39). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates a dense 9 rake-style prongs layout customizable across two distinct tip profiles: Blunt tips for safe tissue manipulation minimizing trauma, or Sharp tips for highly reliable and firm tissue engagement in dense fibrous planes. The ergonomic flat handle ensures fatigue-free control for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Miller Rake Retractor to withstand repetitive hospital sterilization cycles, providing cosmetic, plastic, and reconstructive surgery centers worldwide with exceptional durability, corrosion resistance, and complete autoclave dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)