



RETRACTOR MAKER

PRODUCT BROCHURE OF

Williams Discectomy Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Williams Discectomy Retractor

SKU RM-WDR-1

\$94 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
50 mm × 20 mm (Left or Right)	RM-WDR-1-1	\$94
70 mm × 20 mm (Left or Right)	RM-WDR-1-2	\$99

PRODUCT DESCRIPTION

The Williams Discectomy Retractor is a precision spinal surgical retractor engineered specifically to establish and maintain an unhindered micro-corridor to the interlaminar space during lumbar discectomy procedures. Designed with an angled bayonet-style blade shank, this retractor keeps the assistant's hand and handle completely outside the microscopic line of sight, allowing uninterrupted visualization of the disc space.

Built for targeted paraspinal tissue holding, this Williams Discectomy Retractor model features an overall reach of 18 cm. The instrument is available in two selectable blade sizes: 50 mm × 20 mm (\$94) and an extended 70 mm × 20 mm (\$99) profile, paired with customizable Left-Oriented or Right-Oriented blade directions to support specific unilateral surgical approaches. Its distal toothed tip seats firmly against the laminar border or facet capsule, preventing accidental tissue slippage during deep canal decompressions.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Williams Discectomy Retractor to strict clinical tolerances, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for spine surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)