



RETRACTOR MAKER

PRODUCT BROCHURE OF

Balfour Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Balfour Retractor

SKU RM-BR-STD-1

\$299 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Balfour Spreading Frame	200 × 155 mm, 200 × 235 mm, or 200 × 265 mm
1 Pair	Lateral Retraction Blades	Integrated/Detachable (60 mm or 90 mm depth)
1	Adjustable Center Blade	47 × 80 mm, 62 × 76 mm, or 80 × 80 mm
1	Protective Storage Package	Individual transit box

PRODUCT DESCRIPTION

The Balfour Retractor is a foundational self-retaining surgical retractor engineered to provide stable, hands-free wound edge separation during general abdominal and pelvic operations. By utilizing a smooth sliding crossbar and positive screw lock mechanism, this instrument maintains bilateral abdominal wall retraction while minimizing assistant fatigue.

Designed to accommodate a wide variety of surgical incisions and patient physical profiles, the Balfour Retractor is manufactured in three standardized clinical configurations: a 200 × 155 mm frame (with 47 × 80 mm center blade and 60 mm deep lateral blades), a 200 × 235 mm frame (with 62 × 76 mm center blade and 60 mm deep lateral blades), and a large 200 × 265 mm frame (with 80 × 80 mm center blade and 90 mm deep lateral blades). Each option features an adjustable third-arm center blade to displace the bladder or visceral contents atraumatically.

Manufactured from high-tensile, Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin finish, Retractor Maker supplies the Balfour Retractor non-sterile and fully autoclavable inside individual protective packaging, delivering rugged durability and reliable clinical performance to surgical suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)