



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cottle Alar Nasal Retractor / Hook

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cottle Alar Nasal Retractor / Hook

SKU RM-CANR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Ball on Left Side — 14.0 cm length, 10 mm width, left-sided ball prong	RM-CANR-1-1	\$44
Ball on Right Side — 14.0 cm length, 10 mm width, right-sided ball prong	RM-CANR-1-2	\$44

PRODUCT DESCRIPTION

The Cottle Alar Nasal Retractor / Hook is an essential rhinoplasty surgical instrument engineered specifically to provide safe, non-slip retraction of the alar rim, nostril margin, and lower lateral cartilages. Designed with a distinctive hybrid dual-prong architecture, this instrument utilizes one blunt ball-ended prong alongside one sharp prong to deliver optimal mechanical traction without puncturing fragile cutaneous surfaces.

Built for targeted anatomical alignment, this Cottle Alar Nasal Retractor / Hook model features an overall length of 5 ½ inches (14.0 cm) and a 10.0 mm wide working end (\$44). Surgeons can select between two dedicated operational orientations: Ball on Left Side or Ball on Right Side, ensuring that the protective ball prong faces the external skin margin while the sharp prong anchors securely into the fibrous vestibular lining. Its flat ergonomic handle provides refined manual tactile feedback.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cottle Alar Nasal Retractor / Hook to withstand repetitive hospital sterilization cycles, providing plastic and ENT surgical suites worldwide with high corrosion resistance, unyielding precision, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)